

FFP2/KN95 Particle Mask


Product Specifications

Ultra-fine fibre electrostatic fabric and PP sub bonded non-woven cloth are used to form 4 filtration layers, which can filter harmful substances effectively.

The three-dimensional shape is designed to fit tightly and offer increased breathing comfort.

Fitting instructions











Scope of applications

Suitable for the protection against harmful particles such as dust, PM2. 5 particulate and influenza bacteria droplets.

KN95/FFP2 standard is designed to filter out at least 94% of airborne particles.

Notes

- Wash hands before fitting the mask.
- To ensure good hygiene, avoid touching the inner side with your hands
- Once fitted, check for a tight seal around nose and chin
- Replace the mask if visible dust occurs or breathing resistance increases
- This mask cannot be washed
- Unused masks should be stored in a clean, dust free environment

FFP2 / KN95

